


Understanding Your Hearing Aid

A hearing aid is a highly technological device that can empower you to hear better days and live your best life. This highly sophisticated communication device is programmed by your hearing care professional to meet your specific needs.

The Miracle-EarENERGY™ RIC R is roughly the size of a quarter.




TURN YOUR MIRACLE-EAR APP INTO A REMOTE CONTROL FOR YOUR HEARING AID

The Miracle-Ear App allows you to easily and quickly adjust the functions of your hearing aid so you can hear optimally in the most diverse situations. Almost all our hearing aid devices are compatible with the Miracle-Ear app, making it easier than ever for you to stay connected to the world around you.




AVAILABLE COLORS

Variety of color options to match your style



DIGITAL CHIP

Also known as the brain of the hearing aid, the digital chip amplifies and analyzes sound to provide a high definition digital sound quality.



BLUETOOTH

Easily listen to music, television and phone calls directly to your hearing aids.



MOTION SENSOR

Whether walking or running, the motion sensor automatically responds to deliver the most personalized hearing experience.



DISCREET

Our most popular RIC model sends sound from the hearing aid directly to the ear through a near piece or custom earmold for the utmost comfort.




ACOUSTIC SENSORS

Understands the acoustic environment to deliver the most optimal listening in any situation.



MICROPHONE

Allows you to focus on conversations and improves speech understanding and noise.



CONTROL

Control all your hearing aid volume and settings directly from your hearing aid for your own listening preferences. Hearing aid can also be adjusted from a Miracle-Ear App, remote control or automatically.



RECHARGEABLE

Say good-bye to batteries. Fully rechargeable offers a convenient portable charger with built in power.